BioSystems

REAGENTS & INSTRUMENTS

Manual de mantenimiento Maintenance manual

ES

EN

MANUAL DE MANTENIMIENTO..3

ESPAÑOL

MANTEINANCE MANUAL.....7

ENGLISH

TESE-00011-01-P1
10-2010

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1. Introduction

In order for the imLD microscope to work correctly, we recommend that the following preventive maintenance operations be performed at least once a year.


2. Check of mechanical components

a. Check the adjustment of the plate using the sample positioning control knobs. The plate must remain in a fixed position.

- If it has to be lubricated, proceed as follows:


- X axis:

- Clean the guides (1) of the bearing using a dry cloth. To access the guides completely, move the plate forward and back. Using a fine brush, apply silicone grease to the bearing guides (1). In order to thoroughly distribute the grease, move the plate. Wipe off any excess grease using a cloth. **Do not** put grease on the teeth of the rack and pinion.

- Y axis:


- Repeat the lubrication operation on the guides of the Y bearing (2). **Do not** put grease on the teeth of the rack and pinion.

b. Check the vertical movement of the plate using the macrometric and micrometric control knobs, thereby assuring that the plate remains in the desired position. If it is not correct, adjust the tension of the control knobs using the tension adjustment ring of the focus control knob (3).

- For lubrication, use the same type of grease as in the preceding operations. Clean the guides (4) using a dry cloth. To access them completely, raise and lower the sample holder, and clean the guides both above and below the


sample holder. Use a fine brush to apply grease lightly to the guides (4). Spread the grease by moving and lowering the sample holder. Wipe off any excess grease.


c. Check the threads of the objectives, and check that the revolver assembly moves correctly.

d. Check the movement of the interpupillary adjustment. The eyepieces must remain in the set position. Otherwise, the tightness must be adjusted. To do so, loosen the bolts (5) and remove the cover with the scale (6).

ENGLISH

- Loosen the screws (7).

- Using circlip pliers, tighten or loosen the nut (8) until the desired level of tension is achieved.

- Re-tighten the screws (7).

- Put on the cover with the scale (6).

3. Check of the optical components

a. Check the integrity of the LED of the fluorescence module.

- Remove one of the objectives, and select the empty position of the revolver assembly.

- Turn on the fluorescence module. Put the potentiometer at maximum.

- Put a sheet of white paper on the plate.

- Lower the carriage to the lowest position.

- Observe that the projected light matches the light in the photograph.

- All the rectangles are very well defined, and the corona that surrounds them, although diffused, is clean, meaning that it shows no running or stains.

b. Check the halogen lamp.

- Clean if necessary. Use a cloth moistened with alcohol. Avoid touching it with your fingers. Fingerprints or stains reduce its useful life.

4. Cleaning

Necessary material:

- A mixture of 70% alcohol and 30% ether.
- Lens cleaning tissue.
- Rubber bulb for blowing air.
- Fine, soft brush.
- Cotton tipped.

Procedure

a. Cleaning the body of the microscope.

- Use a cloth moistened with a solution of water and neutral soap. Never use detergents or abrasives.

b. Cleaning of the microscope's external components.

- General rules

Once the external optical components have been removed (lenses, condenser, objectives), proceed as follows:

- Use the fine, soft brush to eliminate dust and other particles. Use the rubber bulb to assure that all dust is removed. To eliminate adhered dust, use the lens cleaning paper by wiping in a circular motion. If there are any fingerprints or remains of grease, use a cotton swab moistened with the alcohol/ether solution. Caution: Do not use the alcohol/ether solution to clean the daylight filter. Instead, use distilled water and dry with a soft cloth.

c. Cleaning the internal lenses.

- If dirt is detected on the internal lenses, they must be taken off and be cleaned with the fine brush and compressed air (or using the rubber bulb). It is not advisable to remove the optical elements, given that they would have to be adjusted again.

List of cleaning accessories supplied by Bioystems

Code	Description
AC15506	Air blower
AC15504	Lens tissue
AC15505	Cotton tipped
AC15503	Lens cleaning brush
AC15507	Complete cleaning kit: - Dust blower - Lens cleaning paper - Cotton swabs - Brush for lenses

ENGLISH

BioSystems

Manufactured by: BioSystems, S.A.

Costa Brava, 30, 08030 Barcelona - Spain Tel: 34-93 311 00 00 Fax: 34-93 346 77 99
e-mail: biosystems@biosystems.es <http://www.biosystems-sa.com>